


VANGUARD REFERENCE PATCHES

(Suitable for CARTO® 3)



SAFE SENSOR CONNECTION


The clip was developed to prevent the unintended detachment of the sensor e.g. due to patient or cable movement.

SPACE SAVING DESIGN

The patches offer greater flexibility and more space during use, e.g. for additional cables and hoses, and especially for older patients with smaller thoraxes.

EXCELLENT ADHESIVE PROPERTIES

The patches provide safe holding properties for any application situation, even after they have been repositioned.



SET CONSISTS OF SIX
VANGUARD REFERENCE PATCHES

ARTICLE OVERVIEW – Patches are available in packages with the following catheters suitable for CARTO® 3

Article group	Article description catheter	REF CE Package	REF OEM Catheter	Tubing set required*
Package Irrigated Ablation Catheter, 3D-mapping	NAVISTAR® THERMOCOOL® 8F D-Type Uni	35233	34H37M	yes
Package Irrigated Ablation Catheter, 3D-mapping	NAVISTAR® THERMOCOOL® NAV 8F D-F Type Bi	35231	36H35M	yes
Package Irrigated Ablation Catheter, 3D-mapping	NAVISTAR® THERMOCOOL® NAV 8F J Type Uni	35237	34HJ7M	yes
Package Irrigated Ablation Catheter, 3D-mapping	NAVISTAR® THERMOCOOL® 8FB-Type Uni	35234	34H17M	yes
Package Irrigated Ablation Catheter, 3D-mapping	NAVISTAR® THERMOCOOL® 8F F-Type Uni	35236	34H57M	yes
Package Irrigated Ablation Catheter, 3D-mapping	NAVISTAR® THERMOCOOL® NAV 8F F-J Type Bi	35232	36H5JM	yes
Package Irrigated Ablation Catheter, 3D-mapping	NAVISTAR® THERMOCOOL® 8F C-Type Uni	35235	34H27M	yes
Package Irrigated SF Ablation Catheter, 3D-mapping	THERMOCOOL® SF NAV 8F F-J Type Bi	35241	BNI35FJH	yes
Package Irrigated SF Ablation Catheter, 3D-mapping	THERMOCOOL® SF Nav 8F D-F Type Bi	35238	BNI35DFH	yes
Package Irrigated SF Ablation Catheter, 3D-mapping	THERMOCOOL® SF NAV 8F F-Type Uni	35244	D131803	yes
Package Irrigated SF Ablation Catheter, 3D-mapping	THERMOCOOL® SF Nav 8F D-Type Uni	35243	D131802	yes
Package Irrigated SF Ablation Catheter, 3D-mapping	THERMOCOOL® SF NAV 8F B-D Type Bi	35239	BNI35BDH	yes
Package Irrigated SF Ablation Catheter, 3D-mapping	THERMOCOOL® SF NAV 8F D-D Type Bi	35242	BNI35DDH	yes
Package Irrigated SF Ablation Catheter, 3D-mapping	THERMOCOOL® SF NAV 8F J-Type Uni	35245	D131804	yes
Package Irrigated SF Ablation Catheter, 3D-mapping	THERMOCOOL® SF NAV 8F B-F Type Bi	35240	BNI35BFH	yes
Package EP-Diagnostic Catheter, variable/fix, Loop, 3D-mapping	LASSO® NAV eco	35149	D134901	no
Package EP-Diagnostic Catheter, variable/fix, Loop, 3D-mapping	LASSO® NAV eco	35150	D134902	no
Package EP-Diagnostic Catheter, variable/fix, Loop, 3D-mapping	LASSO® NAV eco	35151	D134903	no
Package EP-Diagnostic Catheter, variable/fix, Loop, 3D-mapping	LASSO® NAV eco	35152	D134904	no
Package EP-Diagnostic Catheter, variable/fix, Loop, 3D-mapping	LASSO® NAV eco	35153	D134905	no
Package EP-Diagnostic Catheter, variable/fix, Loop, 3D-mapping	LASSO® 2515 NAV eco	35170	D134302	no
Package EP-Diagnostic Catheter, variable/fix, Loop, 3D-mapping	LASSO® 2515 NAV eco	35169	D134301	no

*All catheters are compatible with both tubing sets CF and Sma.

TECHNICAL DATA

OEM REF-No.	Electrode diameter	Cable length	Packaging unit
CREFP6	80 mm	800 mm	6 pieces



COMPLEMENTARY ARTICLES

<p>TUBING SET FOR COOLFLOW™ PUMP</p> <p>Tubing Set CF for irrigated RF-Ablation-Catheters with Luer connector, length: 280 cm Art.no.: 1000001767 OEM REF: 103-050</p>
<p>TUBING SET FOR SMARTABLATE™ PUMP</p> <p>Tubing Set Sma for irrigated RF-Ablation-Catheters with Luer connector, length: 400 cm Art.no.: 2000010029 OEM REF: 103-523</p> <p>The tubing sets are commercial goods of the original manufacturer MedFact.</p>

DECLARATION OF CONFORMITY

 <p>This product bears a CE mark in accordance with the requirements in the EC directive 93/42/EEC</p>

The NAVISTAR®, THERMOCOOL®, LASSO®, Carto®, CoolFlow™ and SmartAblate™ trademarks are trademarks or registered trademarks of Biosense Webster, Inc. or Biosense Webster EMEA, a division of Johnson & Johnson Medical N.V./S.A., respectively.


<p>Manufacturer: VANGUARD AG Landsberger Straße 266 12623 Berlin, Germany</p>		<p>Quality  Made in Germany</p>	<p>Further information regarding this topic and our products available under: www.vanguard-healthcare.com</p>
--	---	---	--




PRODUCT LABEL


VANGUARD Reference Patches
 geeignet für CARTO® 3


VANGUARD Reference Patches
 geeignet für CARTO® 3

- ☑ Reference Patches geeignet für CARTO® 3
- ☑ Reference Patches suitable for CARTO® 3
- ☑ Parches de referencia aptos para CARTO® 3
- ☑ Patches de référence convient pour CARTO® 3
- ☑ Elettrodo di riferimento adatte per CARTO® 3
- ☑ Referentiepatches geschikt voor CARTO® 3
- ☑ Referanseplaste egnet for CARTO® 3
- ☑ Elektrody odniesienia przeznaczone do systemu CARTO® 3
- ☑ Patches de Referência adequados para CARTO® 3
- ☑ Referenspatchar lämpliga för CARTO® 3



Vanguard AG
 Landsberger Str. 266
 12623 Berlin, Germany
 service@vanguard-healthcare.com

 6x  + 




 (01)0408134343522(17)151223(21)1008569018

REF : 34361
 LOT : CH 2018-000
 SN : 1008569018
 ⚡ : 2019-12



VANGUARD AG, Landsberger Str. 266, 12623 Berlin, Germany